Development of a Next-Generation Antimicrobial Wound Dressing- Results Count

Phil Bowler, MPhil, FIBMS Science & Technology, ConvaTec

24th November 2015

Making a Difference in People's Lives



"/TM all trade marks are the property of their respective owners.
© 2015 ConvaTec Inc.

Overview

- The clinical problem
- The solution
- The evidence



The Clinical Problem – Wound Biofilm

Bacteria exist in two different life-styles

Planktonic

- Single cells / free-living
- In solution
- Fast multiplying high metabolism

Biofilm

- Community of cells
- Attached to an surface
- Encased in a slime (exopolymeric substance, EPS)
- Slow metabolism





Bacteria adapt – the biofilm life cycle

Biofilm Maturation

Planktonic (microscopic, invisible, single cells)

(multi-species, slime, may be macroscopic – visible)

Potential for invading

(planktonic) infection

Microcolony formation

(multi-species, microscopic – invisible)

Biofilm in Nature & Disease

- Biofilm accounts for >80% of all microbial infections in the body²
- 1. Costerton et al. Bacterial biofilms in nature and disease. Ann Rev Microbiol 1987; 41:435-64.
- 2. National Institutes of Health, 2002

ConvaTec

Biofilm is more tolerant to antimicrobials than planktonic bacteria

Example : Biocide susceptibility testing of Serratia marcescens

Data Courtesy of :

Dr Andrew McBain, Manchester Pharmacy School, Univ. Manchester

MIC: **M**inimum Inhibitory **C**oncentration MBC: **M**inimum **B**actericidal **C**oncentration MBEC: **M**inimum **B**iofilm **F**radication **C**oncentration

ConvaTec's Wound Biofilm Science

- 2003. The 10⁵ bacterial growth guideline: Reassessing it's clinical relevance in wound healing. OWM.
- 2004. Biofilms and their potential role in wound healing. WOUNDS.
- 2007. Antimicrobial activity of silver-containing dressings on wound microorganisms using an in vitro biofilm model. Int. Wound J.
- 2008. Assessing the effect of an antimicrobial wound dressing on biofilms. Wound Rep Reg.
- 2009. Clinical experience with wound biofilm and management: A case series. OWM.
- 2012. Potential implication of biofilm in chronic wounds: A case series. JWC.
- 2013. Biofilm delays wound healing: A review of the evidence. Burns & Trauma.
- 2014. A clinical algorithm for wound biofilm identification. JWC.
- 2014. The visualization of biofilms in chronic diabetic foot wounds using routine diagnostic microscopy methods. J Diabetes Res.
- 2014. An in vitro test of the efficacy of an anti-biofilm wound dressing. Int J Pharmaceutics
- 2014. Clinician perceptions of wound biofilm. Int. Wound J.
- 2014. Impact of a novel, anti-microbial dressing on *in vivo*, *Pseudomonas aeruginosa* wound biofilm: Quantitative comparative analysis using a rabbit ear model. *Wound Rep. Reg.*
- 2015. Clinical Biofilms: A Challenging Frontier in Wound Care. Adv Wound Care.
- 2015. Safety and Performance Evaluation of a Next Generation Antimicrobial Dressing in Patients with Chronic Venous Leg Ulcers. Int Wound J.
- 2015. A real-life clinical evaluation of a next-generation antimicrobial dressing on acute and chronic wounds. *J. Wound Care*. (JWC Award winner 2015).
- 2015. Design of a next generation antimicrobial Hydrofiber dressing to combat wound biofilm (submitted to BMJ Innovations).
- 2015. A real-life clinical evaluation of a next-generation antimicrobial wound dressing in the United Kingdom. Submitted to JWC.
- 2015. Clinical safety and effectiveness of a next-generation antimicrobial wound dressing designed to manage exudate, infection & biofilm (submitted to *IWJ*).
- 2015. Management of diabetic foot ulcers using a next generation antimicrobial dressing: A case series (submitted to J Diabetes Research).

Biofilm Delays Wound Healing

- "...biofilms are the principle cause of wound chronicity "1
- "...biofilms represent a fourth major pillar of chronic wound pathogenesis"²
- Biofilm creates a sustained low-grade & ineffective inflammatory response²
- Biofilm impairs epithelial migration & granulation tissue formation²
- Individual bacterial species possess distinct levels of biofilm virulence³
- Multi-species biofilm delays healing more than single-species biofilm⁴
- Bacteria from patients with persistent infections are positive for biofilm formation⁵
- At least 60% of chronic wounds contain biofilm^{6, 7}
- Biofilm delays wound healing: a review of the evidence⁸

^{1.} Wolcott et al. Chronic wounds and the medical biofilm paradigm. J Wound Care 2010;19:45-53

^{2.} Gurjala et al. Development of a novel, highly quantitative in vivo model for the study of biofilm-impaired cutaneous wound healing. WRR 2011;19:400-10

^{3.} Seth et al. Quantitative comparison and analysis of species-specific wound biofilm virulence using an in vivo, rabbit ear model. J Am Coll Surg 2012; 215:388-99

^{4.} Seth et al. Comparative analysis of single-species and polybacterial wound biofilms using a quantitative, in vivo, rabbit ear model. PLoS ONE 2012;7:e42897

^{5.} Sanchez et al. Biofilm formation by clinical isolates and the implications in chronic infections. BMC Infect Dis 2013;13:47

^{6.} James et al. Biofilms in chronic wounds. Wound Repair Regen 2008;16:37-44

^{7.} Kirketerp-Møller et al. Distribution, organization, and ecology of bacteria in chronic wounds. J Clin Microbiol 2008;46:2712-22

^{8.} Metcalf & Bowler. Biofilm delays wound healing: a review of the evidence. Burns Trauma 2013; 1: 5-12.

Is Wound Biofilm Visible?

Suspected biofilm over friable granulation tissue in leg ulcer

Thick opaque biofilm? Unresponsive to antibiotics

Suspected biofilm in non-healing surgical wound

Suspected biofilm in non-healing surgical wound

Visual Signs

Metcalf DG, Bowler PG, Hurlow J. A clinical algorithm for wound biofilm identification J Wound Care 2014; 23: 137-142.

Optimum wound biofilm management today?

• Cleansing

Debridement

 Antimicrobial Agents/Dressings (anti-biofilm dressings?)

Considerations for an Anti-biofilm Wound Dressing

- Breakdown & disrupt biofilm
- **Expose** bacteria for more **effective killing** by antimicrobial agents such as ionic silver
- Prevent biofilm reformation in the wound

- Same base product as AQUACEL[™] Ag / Extra dressings,
 i.e. physical attributes, 1.2% ionic silver
- Two additional components:
 - 1. A *metal chelator* (protects silver and destabilises biofilm)
 - 2. A *surfactant* (reduces surface tension and loosens biofilm surface)
- Synergistic enhancement between ionic silver, surfactant
 & metal chelator* (IP protected)
- * Said et al. Int J Pharm, 2014

The Evidence:

AQUACEL[™] Ag+ Extra[™] Dressing vs Standard Silver Dressings

The Evidence: Clinical evaluation of AQUACEL[™] Ag+ dressing

AOUACEL'Ag Evaluación de producto AQUACEL Ag+ Introduction Detalles de la enfermera y el paciente Paciente (identificador no personal): 469759 Edad: 11 Sexp: F w M Enfermero/a: <u>ALEXANDER businer</u> AMBONER Institución: <u>C.S.ALAMO</u> E-mail: <u>ALEXANDES, byb.in/REMERADIO</u>A-O Central: <u>ALEXANDES, byb.in/REMERADIO</u>A-O Central: <u>ALEXANDES, byb.in/REMERADIO</u>A-O Central: <u>ALEXANDER</u> Julianes Información personal Fecha de inicio de la evaluación: 4-3-2013 te de CanvaTec, su privacidad es importante para nosotr Historia clinica relevante, condiciones concomitantes, test de diagnóstico o tratan nous control repaints, control to the site of the second s iticos.) Inicia de tratamiento en otra Comunidad Autónoma mediante <u>ατοριτίζου μέτου το ποτοποιτο στο τοπο τοποιτοστά πο στοποιτο το τοποιτο το ποραστά το τοποιτο το τοποιτο το τοποιτο το ποραστά το τοποιτο το τοποιτο το τοποιτο το τοποιτο το τοποιτο το τοποιτο το ποι το </u> Úlcera de pie diabético Pierra 2 Pielo Sacro D Trocânter D Maléolo D Otras: changed twice a week Results: Longitud mayor 2'5 cm x Longitud menor 2 cm x Profundidad 04 cm 3-6 meses (más de 6 meses 🗆 Estancada :: DeterioroX ción (marque todos los que presente) CalorX Olor D Exudado purule Tejido de granulación friable J Decoloración del tejido de granulación c wor, incluse al porcentaio de telido del lecho de la berida) Granulación 60% Sospecha de biofilm % Moderado 🗶 Alto --

Case Study 3: Venous leg ulcer Author: Agneta Bergsten RN, Ola Arkbro RN Södersiukhuset Stockholm, Sweden

WP-SH is an 86 year old female that got her first venous leg ulcer 1999. The wound has threads in all years of prime and prime in the terminological of the control interaction of the con

The earlier protocol of care has included dressing changes twice a week with clensing and debridement. Iodosorb has been used locally in the wound to easy the debridement of the necroic tissue. The wound bed has then been dressed with antimicribial silver dressings and a superabsorbent cover dressing. Compression therapy with a compression wrap has been used all the time.

Methodology: On 8/2/2013 a new dressing evaluation was initiated. The wound was washed, debrided and the neorotic tissue was curited before application of ADUACEL Age dressing as primary dressing. Lower dressing and finally wrapped the superstant and a source dressing confined to be whose briting a use it. 010913 After he first dressing change the wound looks cleaner and the neorotic part started to lose. The level of granulation tissue increase and a bridge of epithelialisation can be seen and split the wound in two parts. After 0 weeks treatment the wound area decrease fron 5 x 15 cm to one smaller wound at 3x4 cm and one at 3 x 8. The wound is also covered with granulation tissue and the patient has very little pain in the wound

nare AQUACEL Ag+ progression towards healing is great. The dressing changes frequency has decreased from every third day in the beginning to every sixth day in the end of the test. The patie and her relatives are very happy to see the improvement of the wound during the freatment

Tipo de lesión:

Localización

Tamaño de la Jesión:

Duración de la lesión

Dolor X Biofilm 🕁 Eritema 🗙 Edema 🗙

Necrático 40%

Úlcera vascular🗙

0-3 meses

Mejoría 🖂

Esfacelado ____%

Sana n

Bajo 🖂

Macerada / húmeda🗙

Seca / eccematosa

Walker et al. A real-life clinical evaluation of a next-generation antimicrobial dressing on acute and chronic wounds. J. Wound Care. (JWC Award winner 2015).

Baseline Information

Clinical Evaluation of a Next-Generation Antimicrobial Dressing on Acute and Chronic Wounds

Percentage decrease (
) or increase (
) for each wound in the evaluation

- Average wound closure in 113 wounds was 73%
- 17% of wounds healed completely
- 63% of wounds achieved at least 75% closure
- Average treatment period was 4.1 weeks
- Wounds that increased in size were associated with aggressive debridement of devitalised tissue

-60

-100

-80

• Dressing did not contain strengthening yarn or have the additional absorptive capacity of AQUACEL[™] Ag+ Extra[™] dressing Walker et al. A real-life clinical evaluation of a next-generation antimicrobial dressing on acute and chronic wounds. JWC. 2015; 1.

-40

Clinical case study examples

<u>Lithuania</u>

Post-traumatic leg wound; Purulent exudate, odor, poor granulation tissue, suspect biofilm, despite antibiotics; AQUACEL[™] Ag+: infection resolves, healed in 7 weeks

Netherlands, Alita Jaspar Stalled traumatic leg wound; Peri-wound maceration, malodorous exudate, despite silver sulfadiazine; AQUACEL[™] Ag+: infection resolved, great improvement

<u>Portugal, Dr Vitor Santos</u> Chronic DFU of 6+ months; Odor, exudate, slough, suspect biofilm.

AQUACEL™ Ag+: peri-wound skin improved, wound bed cleared, healed in 5 weeks

Day 45

* Dressing did not contain strengthening yarn or have the additional absorptive capacity of AQUACEL[™] Ag+ Extra[™] dressing

Concluding remarks

- Biofilm exacerbates chronicity & impedes healing in most infectious diseases
- If a wound (acute or chronic) is not healing and is unresponsive to antimicrobial therapies, biofilm is likely to be implicated
- A comprehensive protocol-of-care is likely to be most effective in overcoming wound biofilm (debridement, cleansing, antimicrobial)
- The first wound dressing designed to combat biofilm (within a protocol-of-care) is now available to patients in Europe and has been associated with outstanding clinical outcomes

